

View Online at <https://aerobasegroup.com/nsn/6550-01-031-6435>

Usage Design:

For detecting oxidase producing microorganisms

Composition Ingredients:

Para-amino dimethylaniline

Primary Container Type:

Tube, dispensing

Primary Container Content Quantity:

30.000 disk

Total Primary Container Quantity:

6

Product Name:

Differentiation disk microorganism para-amino dimethylaniline300

Special Features:

6 vials of 30 disks per vial

Material:

Paper

Supplementary Features:

D-day item for dod

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A10000