

NSN 6550-01-119-7249

And Urobilinogen In Urine Protein Ph Leukocytes Ketone Glucose Blood Bilirubin Test Strips And Color Chart

View Online at <https://aerobasegroup.com/nsn/6550-01-119-7249>

Features Provided:

Desiccant agent

Primary Container Type:

Any acceptable

Product Name:

Test strips&color bilirubin bld glucose ketone leukocytes ph

Special Features:

Reading from top to bottom, the nine reagent areas shall be located in the following order: leukocytes, nitrite, ph, protein, glucose, ketone, urobilinogen, bilirubin, and blood; final results shall be obtained in 60 seconds or less except for leukocytes which can be read at 120 seconds; for semiquantitative determinations

Material:

Plastic

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A10000