

# NSN 6550-01-119-7249

And Urobilinogen In Urine Protein Ph Leukocytes Ketone Glucose Blood Bilirubin Test Strips And Color Chart

View Online at <https://aerobasegroup.com/nsn/6550-01-119-7249>

**Features Provided:**

Desiccant agent

**Primary Container Type:**

Any acceptable

**Product Name:**

Test strips&color bilirubin bld glucose ketone leukocytes ph

**Special Features:**

Reading from top to bottom, the nine reagent areas shall be located in the following order: leukocytes, nitrite, ph, protein, glucose, ketone, urobilinogen, bilirubin, and blood; final results shall be obtained in 60 seconds or less except for leukocytes which can be read at 120 seconds; for semiquantitative determinations

**Material:**

Plastic

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A10000