

View Online at <https://aerobasegroup.com/nsn/6550-01-200-5769>

Usage Design:

Human chorionic gonadotropin antibody determination

Test Capability Quantity:

100

Component Breakdown Within Each Test Kit:

One 10 ml bt antibody conjugate; one 20 ml bt zero diluent calibrator; one 4 ml bt positive reference alibrator; one 4 ml bt positive control; one 40 ml bt buffer; 100 beads

Supplementary Features:

Medsep item for dod, hybritech inc

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0