

View Online at <https://aerobasegroup.com/nsn/6550-01-200-5769>

**Usage Design:**

Human chorionic gonadotropin antibody determination

**Test Capability Quantity:**

100

**Component Breakdown Within Each Test Kit:**

One 10 ml bt antibody conjugate; one 20 ml bt zero diluent calibrator; one 4 ml bt positive reference alibrator; one 4 ml bt positive control; one 40 ml bt buffer; 100 beads

**Supplementary Features:**

Medsep item for dod, hybritech inc

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0