

View Online at <https://aerobasegroup.com/nsn/6550-01-208-2381>

Usage Design:

An ultra violet method for alcohol determination

Test Capability Quantity:

100

Test Media:

Blood or serum or plasma

Component Breakdown Within Each Test Kit:

One hundred 3 ml bt lyophilized nad-adh enzyme; two 5 ml bt ethynol standard; three 110 ml bt glycine buffer

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0