

View Online at <https://aerobasegroup.com/nsn/6550-01-230-3097>

Usage Design:

For use w/the tdx system

Test Capability Quantity:

100

Product Name:

Test kit quinidine detection sufficient for 100 tests

Test Media:

Serum

Component Breakdown Within Each Test Kit:

One 3 ml bottle serum; one 3 ml bottle serum; one 3 ml bottle tracer; one 3 ml bottle pretreatment; 100 sample cartridges; 100 cuvettes

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0