

View Online at <https://aerobasegroup.com/nsn/6550-01-233-7561>

**Usage Design:**

For use w/tdx fluorescence polarization analyzer

**Test Capability Quantity:**

100

**Protective Storage Feature:**

Requires refrigeration between 2 degrees to 8 degrees c

**Product Name:**

Test kit alcohol determination sufficient for 100 determinations

**Test Media:**

Blood or plasma or serum or urine

**Component Breakdown Within Each Test Kit:**

1 test pack; 100 sample cartridges; 100 cuvettes

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0