

View Online at <https://aerobasegroup.com/nsn/6550-01-233-7582>

Usage Design:

For quantitative procainamide determination

Test Capability Quantity:

100

Protective Storage Feature:

Requires refrigeration between 2 degrees to 8 degrees c

Test Media:

Serum or plasma

Component Breakdown Within Each Test Kit:

One 3 ml bt pretreatment solution; one 3 ml bt procainamide antiseru (rabbit); one 3 ml bt procainamide fluorescein tracer) ; 100 sample cartridges; 100 cuvettes

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0