

View Online at <https://aerobasegroup.com/nsn/6550-01-243-2658>

Usage Design:

Rapid id of pathogenic neisseria species: n. Gonorrhoeae, n. Meningitidis, and n. Lactamica

Test Capability Quantity:

25

Product Name:

Test kit neisseria species id 25 tests

Component Breakdown Within Each Test Kit:

Single test tube procedure utilizing a of inoculum 5-10- colonies with set-up and completion within 40 minutes

Supplementary Features:

D-day item for dod

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0