

View Online at <https://aerobasegroup.com/nsn/6550-01-249-8888>

Usage Design:

To detect chlamydia trachomatis in female urogenital swab specimens

Test Capability Quantity:

100

Product Name:

Test kit chlamydia trachomatis detection sufficient for 100tests

Component Breakdown Within Each Test Kit:

100 urogenital swab; 100 treated beads; one 20 ml vial conjugate; one 4 ml vial positive control; one 15 ml vial negative control; one 20 ml vial chlamydia antibody; one 100 ml vial specimen dilution buffer; one bt of 10 tablets of opd (0-phenylenediamine); one 55 ml bt of opd diluent; one 55 ml bt of opd diluent; 1 vial sulfuric acid; 2 racks of 60 tubes; 2 trays of 60 wells

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0