## NSN 6550-01-257-6391

**Physical Form:** 

Modified Sternheimer-malbin Stain - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6550-01-257-6391

Liquid
Usage Design:
For use with the kova system which facilitates differentiation of cellular elements in urine and enhances the visualization of casts in urine
Primary Container Type:
Bottle
Primary Container Content Quantity:
25.000 milliliters
Total Primary Container Quantity:
3
Product Name:
Sternheimer-malbin stain modified 25ml 3s
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No
Fiig:
A10000