

View Online at <https://aerobasegroup.com/nsn/6550-01-267-1437>

Usage Design:

Determines lupus

Test Media:

Serum

Component Breakdown Within Each Test Kit:

Aone 2.5 ml vi rythrotex reagent; one 0.5ml vi positive human control serum; one 0.5ml vi negative human control serum; 50 applicator sticks; 50 capillary tubes; 1 glass slide w/6 separate rings

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0