

View Online at <https://aerobasegroup.com/nsn/6550-01-325-0105>

Usage Design:

Use in qualitative analysis of benzodiazepine in urine

Test Capability Quantity:

40

Protective Storage Feature:

Requires refrigeration between 2 degrees to 8 degrees c

Component Breakdown Within Each Test Kit:

80 prefilled vials of sheep antibodies; 44 result cards

Supplementary Features:

Medsep item for dod, syva company

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0