

View Online at <https://aerobasegroup.com/nsn/6550-01-349-4991>

Component Breakdown Within Each Test Kit:

Naphthol as-bi phosphoric acid, 10 ml; fast garnet gpc base solution, 10 ml; tartrate solution, 10 ml; acetate solution, 50 ml; sodium nitrite solution, 10 ml; citrate solution, 50 ml; hematoxylin solution, gill no.3, 50 ml;

Special Features:

Disposition after initial use: nonreusable;

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0