

View Online at <https://aerobasegroup.com/nsn/6550-01-354-9460>

**Physical Form:**

Liquid

**Usage Design:**

Harris formula used for staining the nuclei of all surgical specimens

**Primary Container Type:**

Bottle, plastic

**Primary Container Content Quantity:**

1.000 gallons

**Product Name:**

Hematoxylin stain liquid f/staining nuclei of sur specimen 1 gal

**Supplementary Features:**

Medsep item for dod, surgipath medical industries inc

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A10000