

View Online at <https://aerobasegroup.com/nsn/6550-01-358-3821>

Component Quantity:

3

Usage Design:

A quantitative analysis of amphetamines in human urine; used on syva ets analyzers, monoclonal amphetamine/methamphetamine assay

Test Capability Quantity:

230

Test Media:

Urine

Component Breakdown Within Each Test Kit:

One 6 ml vial reagent a; one 6 ml vial reagent b; one 13.3 ml vial buffer concentrate

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0