Prussian Blue Staining Kit - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6550-01-359-3449

Component Quantity:

3

Nonsupply Items And Quantities:

One 250 ml bt potassium ferrocyanide solution; one 250 ml bt hydrochloric acid solution; one 30 ml bt paraosaniline solution

Special Features:

For hematologic stains; may be useful in study of iron deficiency anemia hemachromatosis and hemosiderosis

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800