

View Online at <https://aerobasegroup.com/nsn/6550-01-360-5724>

Component Quantity:

3

Nonsupply Items And Quantities:

250 ml bt ziehl-neelsen carbol fuschin 1; 250 ml bt of 1 pct. Sulfuric acid 1; 250 ml bt loeffler's methylene blue 1

Usage Design:

Used to stain acid-fast bacteria in tissue sections

Product Name:

Fite stain kt used to stain acid-fast bacteria in tissue section

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800