

View Online at <https://aerobasegroup.com/nsn/6550-01-362-5330>

Usage Design:

For use w/bio-rad diamat clinical chemistry analyzer

Test Capability Quantity:

1000

Component Breakdown Within Each Test Kit:

4 analytical column pre-filters; 2 bt elution buffer no. 1; 1 bt elution buffer no. 2; 2 bt elution buffer no. 2; 2 bt elution buffer no. 3; 1 bt wash solution; ten 10ml reference standard; diluent; hamolysis reagent; 1000 sample cups

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0