

View Online at <https://aerobasegroup.com/nsn/6550-01-369-0514>

Usage Design:

For cardiovascular surgery, catheterization laboratories (pt cas), and vascular surgery (automated coagulation timer cartridges)

Total Primary Container Quantity:

50

Protective Storage Feature:

Requires refrigeration between 2 degrees to 8 degrees c and expiration dated

Test Media:

Blood

Exterior Container Quantity:

100

Special Features:

Requires a total of 800 ul fresh whole blood for duplicate determinations and are stable for 24 weeks when stored at 2-25 celcius degrees;

50 dual channel cartridges per box; 2 boxes per pg

Intermediate Package Quantity:

2

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0