

View Online at <https://aerobasegroup.com/nsn/6550-01-465-5104>

Nonsupply Items And Quantities:

High glucose control solution (blue cap) 3 ml volume, 1; low glucose control solution (white cap) 3 ml volume, 1

Sterility:

Sterile

Disposition After Initial Use:

Disposable

Special Features:

Medisense precision glucose control solutions; for in vitro diagnostic use; designed for use with precision blood glucose system, medisense 2, and comparion 2 blood glucose testing systems; reagent composition: high glucose control solution, 0.22% w/v glucose, 99.70% w/v non-reactive ingredients; low glucose control solution, 0.04% w/v glucose, 99.46% w/v non-reactive ingredients

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800