NSN 6550-01-465-8793

Differentiation Kit - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6550-01-465-8793

Р	ro	d	u	ct	Ν	ar	ne:	•

Kit differentiation use f/bone marrow smears fna's & parasites

Special Features:

Contains 500 ml bottle fixative reagent; 500 ml bottle eosin reag ent; 500 ml bottle methylene blue reagent; used for bone marrows, smears, fna's and parasites such as helicobacter pylori. A quick alternative to traditional wright's and giemsa staining procedures.

Shelf Life:

N/a

Unit Of Measure:

__

Demilitarization:

No