

View Online at <https://aerobasegroup.com/nsn/6550-01-479-1785>

Protective Storage Feature:

Requires refrigeration between 2 degrees to 8 degrees c

Usage Design:

For determining relative sensitivity of microorganisms for susceptibility to penicillin g

Test Capability Quantity:

250

Composition Ingredients:

Penicillin g

Primary Container Type:

Tube, dispensing

Primary Container Content Quantity:

50.000 disk

Total Primary Container Quantity:

5

Product Name:

Sensitivity discs diagnostic penicillin g 250s

Active Ingredients Per Unit:

10.0000 units per disk

Special Features:

Each cartridge contains 50 discs; there are 5 cartridges per package; total of 250 discs

Material:

Paper

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A10000