

View Online at <https://aerobasegroup.com/nsn/6550-01-479-1785>

**Protective Storage Feature:**

Requires refrigeration between 2 degrees to 8 degrees c

**Usage Design:**

For determining relative sensitivity of microorganisms for susceptibility to penicillin g

**Test Capability Quantity:**

250

**Composition Ingredients:**

Penicillin g

**Primary Container Type:**

Tube, dispensing

**Primary Container Content Quantity:**

50.000 disk

**Total Primary Container Quantity:**

5

**Product Name:**

Sensitivity discs diagnostic penicillin g 250s

**Active Ingredients Per Unit:**

10.0000 units per disk

**Special Features:**

Each cartridge contains 50 discs; there are 5 cartridges per package; total of 250 discs

**Material:**

Paper

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A10000