

View Online at <https://aerobasegroup.com/nsn/6550-01-479-4513>

**Usage Design:**

To determine alcohol concentrations in saliva which correlates with blood alcohol concentrations

**Test Capability Quantity:**

24

**Total Primary Container Quantity:**

1

**Product Name:**

Alco-screen

**Test Media:**

Saliva

**Component Breakdown Within Each Test Kit:**

24 individual foil wrapped plastic strips with reactive tip and a color chart

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0