

View Online at <https://aerobasegroup.com/nsn/6550-01-487-2908>

Component Breakdown Within Each Test Kit:

96 rubella antigen coated microwell strips; 1-50 ml sample diluent bottle; 11 ml anti-human igm conjugated to horseradish peroxidase; 0.3 ml each vial negative calibrator reference standard set; 11 ml peroxide substrate solution a; 11 ml tmb chromogen solution b; 100 ml-1 vial washing buffer; one each well holder

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0