

View Online at <https://aerobasegroup.com/nsn/6550-01-500-5188>

Component Quantity:

51

Product Name:

Testkit, ict filariasis

Special Features:

Kit consists of 25 individually packaged test cards, 25 capillaries and 1 product insert; blood samples collected in a capillary should be used immediately; venous blood can be used for up to 3 days if stored at 2 degrees to 8 degrees c; serum or plasma can be stored for up to 2 weeks at 2 degrees to 8 degrees c

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800