

View Online at <https://aerobasegroup.com/nsn/6550-01-506-7228>

**End Application:**

Intended for rapid qualitative detection of group a streptococcal antigen from throat swabs or confirmation of presumptive group a streptococcal colonies recovered from culture; clia waived

**Usage Design:**

For rapid qualitative detection of group a streptococcal antigen from throat swabs or confirmation of presumptive group a streptococcal colonies recovered from culture

**Test Capability Quantity:**

50

**Product Name:**

Quickvue dipstick strep a test; 50s

**Component Breakdown Within Each Test Kit:**

Testrip, positive control solution, negative control solution, extraction solution a, extraction solution b

**Special Features:**

Lateral flow method; accuracy at 95%, read results at 5 minutes; sensitivity: 92%, specificity: 98%, specimen type: throat swab; positive and negative internal procedure controls; positive and negative liquid external controls; store at room temp; dipstick format; no instrumentation needed; all equipment to perform 50 tests included

**Supplementary Features:**

Clia waived

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0