

View Online at <https://aerobasegroup.com/nsn/6550-01-506-7228>

End Application:

Intended for rapid qualitative detection of group a streptococcal antigen from throat swabs or confirmation of presumptive group a streptococcal colonies recovered from culture; clia waived

Usage Design:

For rapid qualitative detection of group a streptococcal antigen from throat swabs or confirmation of presumptive group a streptococcal colonies recovered from culture

Test Capability Quantity:

50

Product Name:

Quickvue dipstick strep a test; 50s

Component Breakdown Within Each Test Kit:

Testrip, positive control solution, negative control solution, extraction solution a, extraction solution b

Special Features:

Lateral flow method; accuracy at 95%, read results at 5 minutes; sensitivity: 92%, specificity: 98%, specimen type: throat swab; positive and negative internal procedure controls; positive and negative liquid external controls; store at room temp; dipstick format; no instrumentation needed; all equipment to perform 50 tests included

Supplementary Features:

Clia waived

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A510v0