

View Online at <https://aerobasegroup.com/nsn/6550-01-506-8217>

Usage Design:

For rapid detection of group a streptococcal antigen directly from throat swabs and beta hemolytic colonies recovered from culture

Test Capability Quantity:

25

Component Breakdown Within Each Test Kit:

Individually packaged test cassettes rabbit anti-strep a antibody and heat inactivated strep a antigen; extraction reagent a; extraction reagent b; sterile throat swabs; tubes and tips, positive control; negative control; package insert, and procedure card

Special Features:

Rapid test; clearly distinguishable +/- endpoint and two reagent extraction; 5 minute results; o/a accuracy of test is 98%/ sensitivity of 95%, specificity type 86%; built in controls satisfy clia requirements for daily qc; store at room temp; contains positive and negative controls; 25 tests

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0