

View Online at <https://aerobasegroup.com/nsn/6550-01-506-8217>

**Usage Design:**

For rapid detection of group a streptococcal antigen directly from throat swabs and beta hemolytic colonies recovered from culture

**Test Capability Quantity:**

25

**Component Breakdown Within Each Test Kit:**

Individually packaged test cassettes rabbit anti-strep a antibody and heat inactivated strep a antigen; extraction reagent a; extraction reagent b; sterile throat swabs; tubes and tips, positive control; negative control; package insert, and procedure card

**Special Features:**

Rapid test; clearly distinguishable +/- endpoint and two reagent extraction; 5 minute results; o/a accuracy of test is 98%/ sensitivity of 95%, specificity type 86%; built in controls satisfy clia requirements for daily qc; store at room temp; contains positive and negative controls; 25 tests

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A510v0