## NSN 6550-01-542-3783

Blood Grouping Test Reagent Set - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6550-01-542-3783

Mouse igg1 pe

## **Special Features:**

Mouse igg1 pe fluorescence controls are used to estimate the amount of nonspecific binding, as well as binding mediated by interaction with fc receptors, using a bd facs brand flow cytometer; the flow cytometer must be equipped to detect light scatter and the appropriate fluorescence, and be equipped with appropriate analysis software (such as bd cellquest or bd lysys ii software) for data acquisition and analysis; each reagent is supplied in phosphate-buffered saline (pbs) containing gelatin and 0.1 % sodium azide; reagents can be used for immunophenotyping by flow cytometry with a variety of secimen types, including peripheral blood, bone marrow aspirates or biopsies, and other body fluids or tissues; each type of specimen can have different storage conditions and limitations that should be considered prior to collection and analysis

	Life:	

N/a

**Unit Of Measure:** 

--

**Demilitarization:** 

No

Fiig:

A23800