

View Online at <https://aerobasegroup.com/nsn/6550-01-543-4616>

End Application:

For use w/piccolo chemistry analyzer

Product Name:

Panel basic metabolic plus piccolo analyzer 10s

Special Features:

Basic metabolic plus disc panel includes the following tests bun, calcium, chloride, cre, glu, d, ld, mg, na tco2; must be refrigerated at 36-46 deg.F; must not be frozen; expired discs can not be used in the analyzer; storage above 8 deg.C will degrade the panel

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No