

View Online at <https://aerobasegroup.com/nsn/6550-01-573-9614>

Usage Design:

Rapid purification and concentration of high-quality dna from agarose gel

Test Capability Quantity:

50

Product Name:

Zymoclean gel dna recovery kit

Component Breakdown Within Each Test Kit:

50 ml agarose dissolving buffer; 6 ml dna wash buffer (24 ml of ethonol must be added to the dna washer buffer); 50 zymo-spin i columns;

50 collection tubes; 1 instruction manual

Special Features:

Eluted dna suited for use in dna ligation, sequencing, labeling, and pcr

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0