

View Online at <https://aerobasegroup.com/nsn/6550-01-573-9614>

**Usage Design:**

Rapid purification and concentration of high-quality dna from agarose gel

**Test Capability Quantity:**

50

**Product Name:**

Zymoclean gel dna recovery kit

**Component Breakdown Within Each Test Kit:**

50 ml agarose dissolving buffer; 6 ml dna wash buffer (24 ml of ethonol must be added to the dna washer buffer); 50 zymo-spin i columns;

50 collection tubes; 1 instruction manual

**Special Features:**

Eluted dna suited for use in dna ligation, sequencing, labeling, and pcr

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0