

View Online at <https://aerobasegroup.com/nsn/6550-01-574-4762>

Usage Design:

For use in blood cell differential counts or as a general cytology stain

Primary Container Type:

Bottle

Primary Container Content Quantity:

500.000 milliliters

Total Primary Container Quantity:

1

Product Name:

Dipquick fixative refill 500ml

Special Features:

Gives coloration similar to standard wright-giemsa stain; rapid 3-step procedure can be done in 15 seconds in widemouth 180 ml container; does not require refrigeration; 500 ml refill bottle

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A10000