

View Online at <https://aerobasegroup.com/nsn/6550-01-574-4762>

**Usage Design:**

For use in blood cell differential counts or as a general cytology stain

**Primary Container Type:**

Bottle

**Primary Container Content Quantity:**

500.000 milliliters

**Total Primary Container Quantity:**

1

**Product Name:**

Dipquick fixative refill 500ml

**Special Features:**

Gives coloration similar to standard wright-giemsa stain; rapid 3-step procedure can be done in 15 seconds in widemouth 180 ml container; does not require refrigeration; 500 ml refill bottle

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A10000