

View Online at <https://aerobasegroup.com/nsn/6550-01-600-2397>

**Special Features:**

Swab suspected stain and insert into pouch. Close pouch and break left ampoule. Agitate pouch, then break right ampoule and agitate again. If mammal blood is present a blue-green color will develop within a few seconds. Left ampoule is water, acetic acid, and less than 1 percent leucomalachite green. Right ampoule is 3 percent  $H_2O_2$ .

**Functional Description:**

Presumptive identification of mammal blood.

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No