

View Online at <https://aerobasegroup.com/nsn/6550-01-600-2397>

Special Features:

Swab suspected stain and insert into pouch. Close pouch and break left ampoule. Agitate pouch, then break right ampoule and agitate again. If mammal blood is present a blue-green color will develop within a few seconds. Left ampoule is water, acetic acid, and less than 1 percent leucomalachite green. Right ampoule is 3 percent H_2O_2 .

Functional Description:

Presumptive identification of mammal blood.

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No