

View Online at <https://aerobasegroup.com/nsn/6550-01-600-8645>

End Application:

Bd cytometric bead array (cba)

Special Features:

Employs particles with discrete fluorescence intensities to detect soluble analytes at very low concentrations. Contains all of the supporting reagents necessary to perform an assay. Buffers and instrument setup reagents provided in this kit have been optimized for analysis of analytes in tissue culture supernatants, plasma, and serum samples. Provides sufficient reagents for the quantitative analysis of 100 samples and 10 instrument setup procedures.

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No