

View Online at <https://aerobasegroup.com/nsn/6550-01-625-5026>

**Usage Design:**

Detects all shiga toxin producing e.Coli

**Test Capability Quantity:**

30

**Protective Storage Feature:**

Requires refrigeration between 2 degrees to 8 degrees c

**Test Media:**

Feces

**Component Breakdown Within Each Test Kit:**

E. Coli shiga toxins 1 and 2

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A510v0