

View Online at <https://aerobasegroup.com/nsn/6630-01-047-5412>

**General Description:**

Reflectance operation method, used for whole blood glucose determination in conjunction w/dextostix reagent strips, 10 mg per 100 ml to 400 mg per 100 ml range, porm 5 mg to porm 10 mg per 100 ml precision, 115 vac, 50/60 hz, 3 w, 180 mm by 110 mm by 50 mm portable case size, 10000 hr lamp life

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

A23900