

View Online at <https://aerobasegroup.com/nsn/6630-01-217-4156>

Usage Design:

Colorimetric analysis

Color Range:

+0.0/+1.0 ppm phenolic compound and +0.0/+12.0 ppm phenolic compound

Primary Container Type:

Ampul

Primary Container Unit Quantity:

30

Total Primary Container Quantity:

1

Special Features:

Used to determine phenolic compound in wastewaters produced in the manufacture of petroleum products; steel; resins; and organic chemical; used with phenols test kit p/n k8012; 4 amino-antipyrine in alkaline solution; sealed under vacuum

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T149-b