

View Online at <https://aerobasegroup.com/nsn/6630-01-241-9669>

**Product Name:**

Wash sol auto analyzer.130gm/l quartency ammonium compound 4 li

**General Description:**

Liquid; 4 liter; 0.130 gm/l quartency ammonium compound w/surfactant; used to flush blood cell analyzer before placing it in the standby mode, used w/6640-01-256-9697 model 131 baker semi-automated blood cell counter also series 150, 170, 500, 700, 5000, and 7000

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No