NSN 6630-01-241-9669

Automated Analyzer Wash Solution - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6630-01-241-9669

Product Name:

Wash sol auto analyzer.130gm/l quartency ammonium compound 4 li

General Description:

Liquid; 4 liter; 0.130 gm/l quartency ammonium compound w/surfactant; used to flush blood cell analyzer before placing it in the standy mode, used w/6640-01-256-9697 model 131 baker semi-automated blood cell counter also series 150, 170, 500, 700, 5000, and 7000

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No