## NSN 6630-01-320-5627

Clinical Chemistry Analyzer Wash Solution - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6630-01-320-5627

Physical Form:
Liquid
Specific Use:
To wash out sample pathway on ciba corning model 278 blood gas analyzer
Primary Container Type:
Bottle
Primary Container Content Quantity:
750.000 milliliters
Product Name:
Wash sol clinical chemical analyzer aqueous salt sol 750 ml
Special Features:
Aqueous salt solution
Shelf Life:
N/a
Unit Of Measure:
Demilitarization:
No
Fiig:
T329-a