Au Reagent Aspiration Ise Needle - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6630-01-336-9308

## **General Description:**

C/o 1 1/2 in. L cres pipet tip connected to a 2 in. L polytetrafluoroethylene tubing which is connected to a 3 in. L piece of metal; mixes diluted serum, plasma or urine samples in spin cup and then analyzes those diluted samples to determine concentration of sodium, potassium or chhloride ions; used w/il monarch chemistry system

Shelf Life:

N/a

Unit Of Measure:

---

Demilitarization:

No