

View Online at <https://aerobasegroup.com/nsn/6630-01-348-3135>

Physical Form:

Liquid

Specific Use:

Intended in the abl 500 for automatic rinse of the measuring system after each measurment or calibration; built-in enzyme tavlet must be mixed with the rinse solution before use

Primary Container Type:

Any acceptable

Primary Container Content Quantity:

140.000 milliliters

Total Primary Container Quantity:

1

Unit Package Type:

Any acceptable

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

T329-a