

View Online at <https://aerobasegroup.com/nsn/6630-01-348-3135>

**Physical Form:**

Liquid

**Specific Use:**

Intended in the abl 500 for automatic rinse of the measuring system after each measurement or calibration; built-in enzyme tablet must be mixed with the rinse solution before use

**Primary Container Type:**

Any acceptable

**Primary Container Content Quantity:**

140.000 milliliters

**Total Primary Container Quantity:**

1

**Unit Package Type:**

Any acceptable

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T329-a